

Product Catalog 2024

Beauty is an attitude



amiea.®

QUALITY MADE IN GERMANY





PIGMENTS	4
DEVICES	12
POWER UNITS	19
CARTRIDGES	20
TAPER LENGTHS OVERVIEW	25
VYTAL SKIN	26

Evolutionline Pigments

Your most intense eyes and lips.

Perfect choice for natural-looking eyebrows, as the colors were optimized for the finest hairstrokes, soft shading, and powder effects. Additionally, those who desire intense, brilliant lips will appreciate the long-lasting and stable results when using Evolutionline pigments. These pigments are particularly suitable for beginners. It has a minimal risk of pigment migration and crisp healed results.













- Hybrid pigments - organic and inorganic components
- Quicker healing due to inorganic pigments
- Internal correction factor helps avoid red or gray shift
- Intense and long-lasting colors due to little or no Titanium Dioxide content
- Almost no pigment migration
- Can be mixed with other amiea pigments

Vivid, brilliant, and long-lasting results
Perfect for the finest hairstrokes & powder effects
Fewer, shorter treatments & quicker healing
REACH compliant

EYEBROW PIGMENTS – EVOLUTIONLINE

*colors are formulated without titanium dioxide!

warm ☀ neutral ○ cold ❄

 SAND FXB715☀	 SAHARA FXB720❄	 OCHER FXB730☀	 CLAY FXB740○
 SIENA FXB760☀	 SEPIA FXB770○	 MUD FXB780❄	 UMBRA* FXB785○
 LAVA* FXB790❄	 SOIL* FXB795❄		













CORRECTIVE PIGMENTS – EVOLUTIONLINE

 YELLOW* FXYE700☀	 ORANGE* FXOR700☀	 OLIVE* FXOL700❄	 KHAKI FXKH700❄
 DILUTION FXDIL25			

LIP PIGMENTS – EVOLUTIONLINE

*colors are formulated without titanium dioxide!

warm ☀ neutral ○ cold ❄

 STAR CLUSTER FXR725☀	 CORAL ZINNIA FXR715○	 ROSÉ PEONY FXR735❄	 CLASSIC POPPY FXR720☀
 RED POSETTIA* FXR730○	 SWEET PEA FXR765○	 NUDE ROSE FXR780☀	 NUDE ORCHID FXR785❄
 LOVELY FRESIA FXR775❄	 FRESH TULIP FXR760☀	 BRIGHT AMARYLIS* FXR740❄	 BURGUNDY DAHLIA* FXR750❄

EVO 2 Pigments

NEW

Authentic Beauty,
what you see is what you get!









EVO 2 ensures true colors from application to healed results. Our mineral-based hybrid formulation blends organic and inorganic pigments, trust and color intensity, alongside corrective elements, ensuring enduring color stability. Its high pigment concentration provides increased color opacity, saving you the need for extra passes. The new formulation also boasts improved viscosity for perfect implantation.



- Inorganic with a drop of organic stabilizers
- Stable healing without any shift of color
- No need for corrective pigments after healing
- Higher contrast and opacity
- Free of Titanium Dioxide and Azo pigments
- No gray or red healing, non-whitening effect
- Easier pigment implantation, improved flow
- Crisp lines
- REACH compliant

EYEBROW PIGMENTS – EVO 2

warm ☀ neutral ○ cold ❄

 WARM BLOND FXB001 ☀	 BLOND FXB002 ○	 RED BROWN FXB003 ☀	 WARM MEDIUM FXB004 ○ ☀
 MEDIUM BROWN FXB005 ○	 ASH BROWN FXB006 ○ ❄	 WARM DARK FXB007 ○ ☀	 DARK BROWN FXB008 ❄
 EXTRA DARK FXB009 ○ ❄	 INTENSE DARK FXB010 ❄		



Organicline Pigments

Trusted quality for limitless creativity.

Formulated without iron oxides, the Organicline is designed for eyebrows, lips, eyeliner, and other treatments such as areola reconstruction or scar camouflage. The impressive variety of colors allows creative mixing in order to achieve unique results. These pigments are easier to work with due to the use of an innovative solvent that facilitates pigment delivery into the skin. The risk of drying out during the procedure and potential pigment migration is also minimized.














- Organic pigments - no Iron Oxides used
- Easier skin implantation
- No drying out
- Minimal risk of red/gray shift or migration
- Can be mixed with other amiea pigments

Saturated, intense colors for bright results
Flawless results from eyebrows to camouflage
Easier application & strong coverage

EYEBROW PIGMENTS – ORGANICLINE

*colors are formulated without titanium dioxide!








warm ☀ neutral ○ cold ❄

 TOFFEE FXB520 ○	 HAZELNUT FXB540 ○	 ICE COFFEE FXB524 ❄	 NOISETTE FXB614 ☀
 WOOD TREASURE FXA600 ☀	 CHESTNUT FXB660 ☀	 DARK CHOCOLATE FXB560 ○	 PERSIAN HORSE FXB670 ☀
 AMAZON FXB680 ○	 ESPRESSO FXB580 ❄	 GEISHA FXB690 ❄	










EYELINER PIGMENTS – ORGANICLINE

 BLACK BEAUTY* FXBL022 ❄	 DEEP BLACK* FXBL020 ❄
--	--

LIP PIGMENTS – ORGANICLINE

 PAPAYA FXR501 ☀	 BRICK RED FXR590 ○	 NOUGAT FXR510 ☀	 TASTY NUDE FXR515 ☀
 LE ROUGE FXR550 ☀	 CHERRY BLOSSOM FXR530 ○	 FUCHSIA FXR576 ❄	 CRANBERRY FXR555 ○
 SWEET PLUM FXR572 ❄	 ROMANTIQUE FXR585 ○	 DRAGON FRUIT FXR580 ❄	 MERLOT FXR570 ❄

CORRECTIVE PIGMENTS – ORGANICLINE

 ORGANIC CAMO 10 FXC510 ○	 ORGANIC CAMO 40 FXC540 ☀	 ORGANIC AREOLA 30 FXA530 ○	 ORGANIC AREOLA 40 FXA540 ○
 ORGANIC AREOLA 50 FXA550 ○	 ORGANIC AREOLA 97 FXA697 ☀	 CORRECTIVE YELLOW FXYE520 ☀	 CORRECTIVE MANDARIN FXR504 ☀
 DILUTION FXDIL25			

DEEP Pigments

Next-Level organic colors

NEW


DEEP are purely organic formulated pigments that feature six new striking shades, specially tailored for dark skin tones. With less Carbon Black and no Titanium Dioxide, guaranteeing both true PMU color tones and longevity without color shift.

- Special pigment line for darker skin tones
- Stable healing without any shift of color
- Pure organic formulation
- Higher color density
- Free of Titanium Dioxide and with a reduces amount of carbon black
- Long lasting in the skin
- Perfect for shading and hairstrokes



EYEBROW PIGMENTS – DEEP

warm ☀ neutral ○ cold ❄



Scalp Pigments

amiea scalp range is a generation of pigments, especially developed for scalp micropigmentation. The new formulation contains a balanced relation of oil and water – considering the scalp anatomy. Achieve very accurate and natural results with top-of-the-line pigments that deliver a perfectly balanced viscosity and steady ink flow.



Details:

- Excellent combination of organic and inorganic pigments
- Without azo-colorants and aromatic amines
- Low Iron Oxide = safe for MRI procedures
- Mid to long-lasting saturated colors
- Extremely low pigment migration
- Stable color tones and saturation

SCALP PIGMENTS

	NATURAL BLACK FXB806		LIGHT BROWN FXB801		MEDIUM BROWN FXB802		DARK BROWN FXB804
	CORRECTOR FXOR800		CAMOUFLAGE FXC800				

Genius Evolve

Unleash your freedom

The Genius Evolve is a versatile PMU machine specifically designed to meet your needs. With the choice between a stroke of 2.5 or 3.0 mm, wired and wireless working and the support of amiea Genius cartridge system, it provides you with maximum control over color implantation and enables astonishing detailed work.

Genius cartridge system

Compatible with: NEXT Control Unit

2.5 or 3.0 mm

stroke to choose from

Hybrid operation

wired or wireless use

2 rechargeable batteries

up to 4 hours wireless working

Highest comfort

lightweight, minimal vibration



Details:

- Adjustable needle depth
- Operating range: 60 – 150 Hz
- Material: Premium-grade aluminum
- Size: 145 x 25.5 mm (with battery),
125 x 25.5 mm (with cable connector)
- Weight: 103 g (with cable connector)
127 g (with battery)
- Made in Germany

artyst H3 Boundless

The definition of flexibility

The flexible, hybrid machine that breaks all bounds and elevates your artistry. Building on the success of the artyst H1, it takes your permanent makeup to the next level. Offering wired and wireless working, two stroke lengths to choose from, and the universal cartridge system, it is all you need to set your talents free.

UNIVERSAL CARTRIDGES

Universal cartridge system
Compatible with: artyst Power Unit
2.5 or 3.0 mm stroke to choose from

Hybrid operation
wired or wireless use
2 rechargeable batteries
up to 4 hours wireless working
Highest comfort
lightweight, minimal vibration

Details:

- Adjustable needle depth
- Operating range: 60 – 150 Hz
- Material: Premium-grade aluminum
- Size: 140 x 25.5 mm (with battery),
120 x 25.5 mm (with cable connector)
- Weight: 103 g (with cable connector)
127 g (with battery)
- Made in Germany

Available colors:



Genius

Boost safety and precision

Crafted to provide ultimate precision and flexibility, Genius PMU machines are a perfect fusion of power, comfort, and safety. Achieve precise results, whether you prefer subtle shading or intricate hairstrokes, with the three stroke options of Genius machines.

Genius cartridge system
Compatible with: NEXT Control Unit

Best choice for:

1.5 mm stroke

Feel every insertion for gentle results

2.0 mm stroke

Perfect hairstrokes and finest lines

2.5 mm stroke

Ideal & precise shading

Details:

- 1.5, 2.0 or 2.5 mm stroke
- Lightweight: 62 g
- Made in Germany

artyst H2 Powerbabe

Fusion of power and ergonomics

Experience a revolutionary fusion of ergonomics and tattoo power. This powerful yet low-vibration PMU machine, with a 3.5 mm stroke, is a collaboration between amiea and Cheyenne design experts. A perfect PMU tool for artists seeking the highest performance with full flexibility when it comes to choosing needle cartridges.




Universal cartridge system
Compatible with: artyst Power Unit or NEXT Control Unit with adapter
3.5 mm stroke

Best choice for:
Shading & pixel techniques
Fine lines
Crisp hairstrokes
All skin types

Details:

- Adjustable needle depth
- Frequency: 70 - 150 Hz
- Material: Premium-grade aluminum
- Weight: 130 g
- Made in Germany

Available colors:



Supreme

Elegance paired with precision

Supreme is a combination of German engineering and long-time expertise in permanent makeup. Crafted with accuracy and designed for excellence, this machine embodies safety, precision, and an unparalleled PMU experience for both you and your clients.

Supreme cartridge system (click or screw)
Compatible with: NEXT Control Unit
2.0 mm stroke

Best choice for:

Truly perfect for hairstrokes and precise lines

Details:

- Material: Premium-grade aluminum
- Size: 22 x 163 mm
- Weight: 80 g
- Made in Germany



artyst H1 PMU Classic

Versatile machine to fit your needs

This exceptional PMU machine stands for precision, security, and silence, setting a distinguished standard in the market. Feel the inspiration and barely a vibration as you create your magic with finesse.

H1 cartridge system
Compatible with: artyst Power Unit
2.0 mm stroke

Best choice for:
Ultrafine hairstrokes
Shading
Pixel techniques

Details:

- Adjustable needle depth
- Frequency: 70 - 150 Hz
- Material: Premium-grade aluminum
- Size: 136 mm x 21 mm
- Weight: 80 g
- Made in Germany

Available colors:




Symphony II

Light and comfortable

The era of the amiea Symphony continues in a new and improved ergonomic design. The Symphony II performs at a much lower vibration in comparison to the previous edition making this machine very silent and ultra-precise. It is now also compatible with all amiea Supreme needle cartridges, making it suitable for artists of all levels.

Supreme cartridge system (click)
Compatible with: NEXT Control Unit
2.0 mm stroke



Details:

- Improved precision and minimized vibration
- More ergonomic design
- Better grip and highest comfort
- Adjustable needle depth
- Premium grade aluminum housing
- Made in Germany

NEXT Control Unit

The intelligent control unit comes with an LED touch display and timer functionality. Its lightweight and ensures maximum comfort. The NEXT Control Unit can be used with or without the amiea foot switch.

Compatible with all amiea PMU machines (except Sense)
Supports amiea Footswitch



Details:

- Cost-effective, powerful control unit
- Built in treatment-timer-functionality
- Lightweight & modern design
- Frequency: 70 - 150 Hz
- Power input: 15 V / 10 W
- Weight: 350 g
- LED display dimensions: 130 x 50 x 110 mm
- Material: Frame from premium grade aluminum, glass (LED display)
- Made in Germany

Artyst Power Unit

A versatile, universal power supply that is small and lightweight yet robust. It can be placed on metal surfaces due to its powerful magnets. Furthermore, it is very easy to use - just plain and simple.

Compatible with all artyst machines
Optional foot switch



Details:

- Compact design, easy to operate
- 9 easy-to-understand speed levels
- Mountable on metal surfaces
- Frequency: 60 - 150 Hz
- Power input: 15 V / 10 W
- Dimensions: 87 x 76 x 28 mm
- Weight: 150 g
- Made in Germany

amiea Classic Cartridges

Developed specially for our Supreme and Genius machines. Amiea cartridges are the right tool to skyrocket your permanent makeup results. You can be sure that you will find the perfect cartridge for any treatment within our range. Finest surgical steel needles within the highest quality cartridge guarantee better pigment implantation and make sure to keep the skin trauma to a minimum, so your clients will love your work even more.



Anti-slip grip
Precise needle guidance
Strict quality control
Safety Membrane prevents backflow of any fluids

Configuration		Symphony I	Supreme* (click or screw) Symphony II* (click)	Genius	Diameter in mm	Application area
NANO NT	•		√	√	0.25 mm	Eyes, brows, lips, areola, scalp
NANO N2	◦	√	√		0.25 mm	Eyes, brows, lips, areola
NANO N1	◦	√	√	√	0.20 mm	Eyes, brows, lips, areola
1-LINER	•	√	√	√	0.40 mm	Eyes, brows, lips
3-LINER	•••	√	√	√	0.30 mm	Eyes, brows, lips, scalp
1-MICRO	•	√	√	√	0.40 mm	Eyes, lips, scalp
3-MICRO	•••		√	√	0.18 mm	Eyes, brows, lips, areola, scalp, scars
3-OUTLINE	•••	√	√	√	0.25 mm	Eyes, brows, lips
3-POWER	•••		√		0.30 mm	Eyes, lips, areola, scars
3-SLOPE		√	√		0.40 mm	Eyes, brows, lips, areola
4-FLAT	√	√		0.35 mm	Eyes, brows, lips, areola
5-ROUND	•••		√		0.30 mm	Eyes, brows, lips, areola
5-SHADER	•••	√	√	√	0.30 mm	Eyes, brows, lips, areola, scars
5-MAGNUM	•••	√	√		0.30 mm	Eyes, brows, lips, areola, scars
5-SLOPE		√	√		0.40 mm	Eyes, brows, lips, areola, scars
5-POWER	•••		√		0.30 mm	Eyes, brows, lips, areola, scars
7-ROUND	•••	√	√		0.30 mm	Eyes, brows, lips, areola, scars
7-POWER	•••		√		0.30 mm	Eyes, brows, lips, areola, scars
9-MAGNUM	•••••	√	√	√	0.30 mm	Eyes, brows, lips, areola, scars
3-NANO SLOPE			√	√	0.25 mm	Eyes, brows, lips, areola, scars
5-NANO SLOPE			√	√	0.25 mm	Brows, lips
7-NANO ROUND	•••			√	0.25 mm	Brows, lips
3-NANO SPREAD	•••			√	0.25 mm	Skin, areola
5-NANO SPREAD	•••			√	0.25 mm	Skin, areola
7-NANO SPREAD	•••••			√	0.25 mm	Skin, areola

Compatibility – amiea

Universal Cartridges

All you need for perfection Based on the long-term experience of the tattoo cartridge experts, artist cartridges offer the highest quality and comfort to help you level up your craft In our wide range you will sure find a configuration to fit your needs and every procedure. These cartridge have the universal docking system and can be used with artist H2 Powerbabe, artist H3 Boundless and most other PMU-machines.

Support low machine vibration
Safety Membrane prevents backflow of any fluids
Approved by the best PMU-artists in the world
Universal docking system compatible with most PMU-machines
Come in 20-piece boxes





Configuration	Diameter in mm	Taper	Application area
1-NANO N2 0.25 CPT ○	0.25 mm	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola
1-NANO 0.25 LT ○	0.25 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola
1-LINER 0.30 LT ●	0.30 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola
1-LINER 0.30 ST ●	0.30 mm	Short Taper 1.85 mm	Eyes, brows, lips, areola, scars
1-LINER 0.35 LT ●	0.35 mm	Long Taper 4.5 mm	Eyes, brows, lips
1-LINER 0.35 ST ●	0.35 mm	Short Taper 1.85 mm	Eyes, brows, lips, areola, scars
1-MICRO 0.40 ●	0.40 mm	Medium Taper 3.5 mm	Eyes, brows, lips, areola
3-Nano Slope 	0.25 mm	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola, scars
1-LINER 0.40 CPT ●	0.40 mm	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola
3-Liner 0.25 	0.25 mm	Super Long Taper 5.5 mm	Eyes, brows, lips, areola
3-Liner 	0.30 mm	Super Long Taper 5.5 mm	Eyes, brows, lips, areola
5-Flat 	0.35 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola
5-Liner 0.25 	0.25 mm	Super Long Taper 5.5 mm	Eyes, brows, lips, areola, scars
5-Magnum 	0.35 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola, scars
7-Liner 	0.30 mm	Super Long Taper 5.5 mm	Eyes, brows, lips, areola, scars
7-Magnum 	0.30 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola, scars
9-Magnum 	0.35 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola, scars
13-Magnum 	0.35 mm	Long Taper 4.5 mm	Eyes, brows, lips, areola, scars

artyst H1 Cartridges

These PMU-needles are everything you need to achieve perfection with you artyst H1 PMU Classic machine. These cartridges offer the highest quality and comfort to help you level up your craft. Support low machine vibration Safety Membrane prevents backflow of any fluids Exclusively tested by the best PMU-artists in the world.



Support low machine vibration
Safety Membrane prevents backflow of any fluids
Medium soft safety membrane for smooth skinpenetration
Exclusively tested by the best PMU-artists in the world

Configuration	Diameter in mm	Taper	Application area
1-Nano (N2) ◦	0.25	Classic PMU Taper 1.5 mm	
1-Nano Tube •	0.25	Classic PMU Taper 1.5 mm	
3-Nano Slope 	0.25	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola, scars
1-Liner ●	0.40	Classic PMU Taper 1.5 mm	Eyes, brows, lips,areola
3-Outline 	0.25	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola
5-Magnum 	0.30	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola, scars
3-Micro 	0.18	Classic PMU Taper 1.5 mm	Eyes, brows, lips, areola, scars

Capillary & Flow Cartridges

Capillary and Flow cartridges were created with your desires and comfort in mind. The capillary structure within the cartridge tip improves the needle guidance and stores more pigment which means ideal precision, better pigment distribution and longer color flow. You will save your precious time as you will reduce the number of dips and passes needed to achieve perfect results. Quicker color implantation also means less trauma and better skin tolerance which will leave your clients more satisfied.


Higher pigment intake
Up to 10 times longer intervals between dips
Incomparable color flow and precision



Compatibility	Configuration	Universal	Genius	Diameter in mm	Application area
ameia Capillary Flow	1-NANO (2) CAPILLARY FLOW ●		√	0.25 mm	Eyes, brows, lips, areola
ameia Capillary Flow	1-LINER CAPILLARY FLOW ●		√	0.40 mm	Eyes, brows, lips
ameia Capillary Flow	1-MICRO CAPILLARY FLOW ●		√	0.40 mm	Eyes, lips, scalp
artyst Capillary	1-Liner Cap (0,30 LT) ●	√		0.30 mm	Eyes, brows, lips, areola
artyst Capillary	1-Liner Cap (0,30 ST) ●	√		0.30 mm	Eyes, brows, lips, areola, scars
artyst Capillary	1-Liner Cap (0,35 LT) ●	√		0.35 mm	Eyes, brows, lips
artyst Capillary	1-Liner Cap (0,35 ST) ●	√		0.35 mm	Eyes, brows, lips, areola, scars
artyst Capillary	1-Micro Cap (0,40 MT) ●	√		0.40 mm	Eyes, brows, lips, areola
artyst Capillary	1-Nano Cap (0,25 CPT) ◦	√		0.25 mm	Eyes, brows, lips, areola
artyst Capillary	1-NANO CAP (0,25 LT) ◦	√		0.25 mm	Eyes, brows, lips, areola
artyst Capillary	3-Micro Cap (0,18 CPT) ⬢	√		0.18 mm	Eyes, brows, lips, areola, scars

Technical parameters

Taper length and needle diameter



DIAMETER VARIATIONS

0.18 mm
0.25 mm
0.35 mm
0.40 mm



DIAMETER VARIATIONS

0.25 mm
0.30 mm
0.35 mm
0.40 mm



amiea vytal Skin

Cosmetic microneedling on medical level

The first cosmetic microneedling machine on medical level. Our carefully developed microneedling system based on the medical microneedling technology guarantees you maximum safety and accurate skin penetration.

Microneedling cartridge system (screw)
Compatible with: NEXT Control Unit

Details:

- Gentle on skin
- Almost painless experience for the client
- Light and quiet
- Efficient treatment with circular movements
- Low vibration
- Made in Germany

What matters for your client:
No unnecessary skin trauma
Minimal to no downtime
Visible, sustainable, and even results
Excellent comfort and maximum safety



Microneedling safety cartridges

Five different needle configurations enable you to treat your clients according to their individual needs and skin conditions.

Microneedling cartridges available for screw and click machines (amiea vyal Skin, Supreme, Symphony II)

6-needle tilting plate
with 0.3 mm surgical steel needles (screw)

18-needle tilting plate
with 0.3 mm plastic needles (screw and click)

6-needle tilting plate
with 0.7 mm surgical steel needles (screw and click)

6-needle tilting plate
with 1.0 mm surgical steel needles (screw and click)

6-needle tilting plate
with 0.5mm surgical steel needles (screw and click)



- Premium quality surgical steel and plastic needles
- Patented, tiltable needle plate that adapts to skin's surface and limits penetration depth to needle length
- Bio-compatible material
- Patented Safety Membrane to prevent backflow of fluids
- 90° angle to reduce skin injury
- Validated sterilization procedure
- Strictest quality control

amlea®

BERLIN | GERMANY

www.amlea.com
shop.amlea.com

SMA-0147-EN