

PRODUCT DATA SHEET

SYMPHONY I

Lightest and smallest

The Symphony I is our most basic but still professional micropigmentation system. It offers reliable and satisfying results. Many technicians prefer to work with a minimalist device like this, especially when micropigmentation is not the core business of their studio.



APPLICATIONS

Cosmetic micropigmentation for:

- areola / body
- eyebrows, eye-liners, lip-liners and / or filling in lip colour

FEATURES

Precision and power

BENEFITS

- Low epidermal trauma
- Low noise and vibration
- Superior results

Values

- German technology
- R&D and production in Berlin
- Two years warranty

PERFECT RESULTS

Micropigmentation can be used to restore confidence and improve your overall appearance. It can be applied to eyebrows, eyelines and lips.



eyebrows before



eyebrows after



lips before



lips after

TECHNICAL PARAMETERS

Technical parameters

- Type: A-9.1
- Rated voltage: 15 V - (DC)
- Power input: 7 V A max
- Working frequency: Around 120 hits per second
- Drive: Precision DC motor
- Operating mode: Continuous operation

Operating conditions

- Ambient temperature: +50° to +95°F / 10° to 35°C
- Relative humidity: 30% to 75%
- Width x height x depth: 3.94 x 1.77 x 5.51 inches / 10 x 4.5 x 14 cm
- Weight of handpiece: 1.76 oz / 80 g
- Total weight: 0.66 lbs / 10.58 oz / 300 g
- Manufactured in: Germany according to quality standard ISO 9001

Patents

- US 6,505,530; US 6,345,553; US 8,029,527; US 9,504,814
EP 1 495 782; EP 1 618 915
- Additional patents pending