

# amiea

## PRODUCT DATA SHEET

# LINELLE SUPREME

### No. 1 in sales worldwide

The sleek intuitive control unit is digitally controlled. It ensures top precision results which allows you to focus fully on your work.

The amiea Linelle Supreme operates with minimal noise and vibration.



## APPLICATIONS

Cosmetic micropigmentation for:

- areola / body
- eyebrows, eye-liners, lip-liners and / or filling in lip colour

## FEATURES

Adjustable needle speed

## BENEFITS

- High power
- Perfect pain management
- Minimal downtime
- 50 - 150 pps

Needle length adjustable

- Appropriate needle length for each area of application
- Range of depth 0.02 - 0.14 inches / 0.5 - 3.5 mm

Precision and power

- Low epidermal trauma
- Low noise and vibration
- Superior results

Values

- German technology
- R&D and production in Berlin
- Two years warranty

## PERFECT RESULTS

Micropigmentation can be used to restore confidence and improve your overall appearance. It can be applied to eyebrows, eyelines and lips.



eyebrows before



eyebrows after



lips before



lips after

---

## TECHNICAL PARAMETERS

### Technical parameters

- Type: DA30310US
- Rated voltage: 15 V - (DC)
- Power input: 7 V A max
- Working frequency: Selectable from 50 to 150 hits per second
- Drive: Precision DC motor
- Operating mode: Continuous operation

### Operating conditions

- Ambient temperature: +59° to +77°F / 15° - 25°C
- Relative humidity: 30% to 75%
- Width x height x depth: 10.28 x 3.50 x 8.62 inches / 26.1 x 8.9 x 21.9 cm
- Weight of handpiece: 3.53 oz / 80 g
- Total weight: 1.98 lbs / 31.75 oz / 900 g
- Manufactured in: Germany according to quality standard ISO 9001

---

### Patents

- US 6,505,530; US 6,345,553; US 8,029,527; US 9,504,814  
EP 1 495 782; EP 1 618 915; EP 1 958 659; EP 2 462 979
- Additional patents pending