

PRODUCT DATA SHEET

LINELLE II

Most suitable for mobile use

The amiea Linelle II is lightweight and very easy to operate and at the same time providing all the essentials. This basic, reliable device for micropigmentation is perfect for professionals who offer mobile services.



APPLICATIONS

Cosmetic micropigmentation for:

- areola / body
- eyebrows, eye-liners, lip-liners and / or filling in lip colour

FEATURES

Adjustable needle speed

Needle length adjustable

Precision and power

Values

BENEFITS

- High power
- Perfect pain management
- Minimal downtime
- 50 - 140 pps

-
- Appropriate needle length for each area of application
 - Range of depth 0.00 - 0.098 inches / 0.0 - 2.5 mm

-
- Low epidermal trauma
 - Low noise and vibration
 - Superior results

-
- German technology
 - R&D and production in Berlin
 - Two years warranty

PERFECT RESULTS

Micropigmentation can be used to restore confidence and improve your overall appearance. It can be applied to eyebrows, eyelines and lips.



eyebrows before



eyebrows after



lips before



lips after

TECHNICAL PARAMETERS

Technical parameters

- Type: A-9.9-140
- Rated voltage: 15 V - (DC)
- Power input: 7 V A max
- Working frequency: Selectable from 50 to 140 hits per second
- Drive: Precision DC motor
- Operating mode: Continuous operation

Operating conditions

- Ambient temperature: +50° to +77°F / 15° - 25°C
- Relative humidity: 30% to 75%
- Width x height x depth: 3.94 x 1.77 x 5.51 inches / 10.0 x 4.5 x 14.0 cm
- Weight of handpiece: 1.76 oz / 50 g
- Total weight: 0.88 lbs / 14.11 oz / 400 g
- Manufactured in: Germany according to quality standard ISO 9001

Patents

- US 6,505,530; US 6,345,553; US 8,029,527; US 9,504,814
EP 1 495 782; EP 1 618 915
- Additional patents pending